

## ReliZyme™ HFA403

### Highly porous hydrophilic enzyme carriers




#### Product Features

Matrix	polymethacrylate
Appearance	perfect spherical white opaque beads
Functional group	amino epoxide
Oxirane content	min. 30 µmol/g wet
Median pore diameter	40 – 60 nm
True density	> 1.1 g/ml
Water retention	60 – 70 %
Temperature stability range	0 – 10 °C
pH stability range	5 - 8
Recommended storage temperature	4 – 6 °C
Storage time limit	6 months

Available grades	S	M
Particle size range	100 – 300 µm	200 – 500 µm

#### Pressure drop in water at 25°C



#### Remarks

Wear safety glasses and gloves during the resin handling.  
Refer to the Safety Data Sheet for all details.

#### Notice:

All the data and suggestions made herein are based upon our research and are believed to be accurate. However, no guarantee is made or implied since conditions and methods of use of our products are beyond our control. Our products are sold on the conditions that the user will evaluate himself, as well as our formulae and recommendations, to determine their suitability for his own purpose. Also, statements as to the use of our products are not be construed as recommendations for their use in the infringements of any patent.